



Description: Quick Start for configuring and running the Iconics Custom Security Login

OS Requirement: The operating system required is Win XP/2000/NT.

General Requirement:

Iconics Custom Security Login is an application that will link the user with the Iconics security login by way of the SecuGen EyeD OptiMouse mouse. The ISA-Automation toolkit is used to identify and store new and existing users.


Capabilities:

- The Iconics Custom Security Login is an alternative to the standard Security Login. It allows for immediate login to the Iconics security server for existing users
- Using this new tool the security administrator will have the ability to add new users.


Configure Iconics Custom Security Login

To configure the Iconics Custom Security Login, the local machine must have a registered OnlyU.Toolkit along with Genesis32 and MSAccess.


1. The SecuGen EyeD OptiMouse and drivers must be properly installed.
2. Then install ISA-Automation OnlyU application (Administrator setup)




3. After installation you will need to set up the OnlyU simple search example. Go to Start > Program Files > OnlyU > Configuration Tool. The dialog below will show and the user should select the mouse.



4. On the next dialog the user will select the USB device click next and then database connection string must be built as shown below.



5. After the connection has been create and tested successfully, the Simple search behaviors can be set.





6. Check Simple Search and then Setup and for the User Table Name select Tab_Employees and for Biometric column #1 select biodata 1 and for Biometric column #2 select biodata 2
7. All else can be left as default and Finish.
8. Register the OnlyU.Toolkit by going to Start > Program Files > OnlyU > Registration. From there select OnlyU ToolKit.
9. **Employees.mdb must be edited; a "password" column must be added.**
10. Check to make sure that the new Employees.mdb is located in C:\ProgramFiles\ISA\OnlyU\
11. Genesis32 must also be installed on the local machine.

The login utility is ready to be used.

In order create and add new users to the security server along with the fingerprint database:

1. To add new users to the security application you must have administrative permissions for the Iconics security.
2. Input data for the new user in the start window. (The user must be present to input fingerprint data.) Figure 1
3. Click Enroll fingerprints, follow instructions.
4. Now click add user.
5. Iconics Security System Administrator needs to login Figure 2 and the new user password must be confirmed and finally the user can be added. Figure 3

Iconics Security Login - v1.0.0

Firstname:

Lastname:

Age:

Job Title:

Password:

Email:

Biometric:

Figure 1 Start window

Administrator Login

User Name:

Password:

Figure 2 Login dialog

Add User Form

Enter the information of for the new user

Username:

Password:

Confirm Password:

IMPORTANT check if the user is an administrator

☐ Administrator ☐ User Cannot Change Password

Figure 3 New user confirmation